

Wall Dispenser 25/100/250/500/1000 ml


Whenever liquid media need to be dosed simply and precisely, our wall dosing devices are the first choice. All device parts that come into contact with the liquids are made of stainless steel, aluminum and Teflon. If different materials are used, we recommend using a separate dosing device for each material. This avoids cleaning work.

Wolfangel wall dosing devices are all prepared for wall mounting – with just two screws, they are easily attached to the wall.

Lift the pump lever, the medium is sucked from the original container into the glass cylinder. The exact quantity can be read immediately on the scale on the glass cylinder. Push the lever down, the just-dosed medium flows through the outlet into your provided container. Wall dosing devices 1000 ml are pneumatically operated.

Wall Dispenser	25	100	250	500	1000
Volume	25 ml	100 ml	250 ml	500 ml	1000 ml
Operated	manual			pneumatically	
Item no.					
With limit stop	80301	80306	80310	80312	80319
Without limit stop	80300	80305	80308	80311	80316

Subject to technical modifications



Wall Dispenser
stainless steel with limit stop



Wall Dispenser
stainless steel without limit stop